

ALMICRO

HOT AIR OVEN



Specification:

Standard double wall fabrication. Inner chamber made of richly anodized fabricated out of thick mild steel sheet duly finished in stored enamel off-white paint or pastel mat finish colors combinations. The gap between the walls is filled with special grade grass wool for heat insulation. Inner chamber is ribbed for placing shelves at convenient levels. Provided with two or three remove able shelves. Special elegant door closing device has been provided.

Nichorome wire heating elements provided on three sides attain quick and uniform heating in range of 50°C to $250^{\circ} \pm 1^{\circ}\text{C}$ controlled by capillary type thermostat.

L-shaped thermo meter is built in type. Control panel is provided with selector switch of high or low rates of power. Thermostat control knob and indicators for mains & thermostat. Supplied with cord and plug. The equipment is suitable to operate on 220V AC 50 Hz single phase.